



SOOO STABLE
IS **VISALYS®** CORE

CORE BUILD-UP & POST CEMENTATION



KETTENBACHDENTAL

Simply intelligent

VISALYS® CORE – EASY TREATMENT

Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT- TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives



DENTIN-LIKE GRINDABILITY




FREE OF BISPHENOL A



Distinguished by the
experts “The Dental
Advisor”

VISALYS® CORE – STABLE RESULTS


Secure bonding with all relevant adhesives



In-house measurements by Kettenbach, shear adhesive bond, test based on DIN EN ISO 29022.

Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Removal depth when grinding various core build-up materials



Tengku Ibrahim Maulana. Grindability analysis of dental core build-up materials under clinical relevant test procedure. Journal of Medical Materials and Technologies (Thesis Series Edition), Vol 1 No 1, 1-21 (2017)

Visalys® Core can be precisely processed similar to dentin.

VISALYS® CORE – THE RANGE



Visalys® Core

Intro pack

White
Dentin

1 x 9 g (5 mL) automix syringe,
10 brown mixing tips,
5 intraoral tips, 5 endo tips

Normal pack

White
Dentin

2 x 9 g (5 mL) automix syringe,
20 brown mixing tips,
10 intraoral tips, 10 endo tips

Normal pack

White
Dentin

1 x 45 g (25 mL) cartridge,
20 yellow mixing tips,
20 intraoral tips



NOTES

Kettenbach SNC
7, Place de la Gare
57200 Sarreguemines · France
info@kettenbach.com
www.kettenbach-dental.com



022740_1024

KETTENBACHDENTAL
Simply intelligent